

Point-of-Care Testing Solutions

Get the results you need when you need them with CardioChek[®] Plus and A1CNow^{®+}





Ø



ptsdiagnostics.com

CardioChek® Plus


Key Benefits

- Accurate results in as little as 90 seconds
- Portable
- Test strip refrigeration not required
- CLIA-waived


Accuracy

- Meets NCEP guidelines for accuracy and precision
- CRMLN-certified for CHOL and HDL tests
- FDA-cleared


Cholesterol Linear Regression



Triglycerides Linear Regression




HDL Cholesterol Linear Regression



Roche Cobas Integra 400 Plus HDL Cholesterol (mg/dL) Average Bias = -4.67%

Glucose Linear Regression



Product Specifications

Measures

- Total Cholesterol
- HDL Cholesterol
- Triglycerides
- Glucose

Calculates

- LDL Cholesterol
- TC/HDL Ratio
- LDL/HDL Ratio
- Non-HDL Cholesterol

Time to Result

As little as 90 seconds

Sample Volume

- Lipid Panel Fingerstick (40 μL)
- eGLU Fingerstick (1.1 μL)

Storage and Handling

Refrigeration is not required. But product may be refrigerated at room temperature 68-86 °F (20-30 °C) or refrigerated 35-46 °F (2-8 °C). Good until stated expiration date.

Power Source

4 AA batteries

Test Range

Lipid Panel:

Total Cholesterol	100-400 mg/dL
HDL Cholesterol	20-120 mg/dL
Triglycerides	50-500 mg/dL
eGLU:	40-600 mg/dL

Methodology

Employs both light reflectance and electrochemical biosensor technology to measure an enzymatic chemical reaction.



A1CNow^{®+}

Key Benefits

- Accurate results in 5 minutes
- Portable
- CLIA-waived

Accuracy


The A1CNow⁺ system has demonstrated accuracy equivalent to clinical laboratory analyzers.

CE-marked

- NGSP-certified
 FDA-cleared
- IFCC-traceable



HbA1c Linear Regression



Accuracy Average Bias =

-0.03 %A1C

Bland-Altman (95% confidence) = -0.6 to +0.5 %A1C

Precision (14 samples / 10 replicates each) Average %CV = 3.34%

Lab Analyzers Used in Evaluations

Beckman Coulter AU5400, Roche Cobas[®] Integra 400 plus, Roche Cobas c513, Roche Cobas Integra 800, Siemens Dimension[®], Tosoh G8



Product Specifications

Measures %A1C

Time to Result 5 minutes

Sample Volume Fingerstick (5 µL)

Storage and Handling

Store refrigerated 36-46 °F (2-8 °C) until expiration date. Stored up to 4 months at room temperature 64-82 °F (18-28 °C).

Power Source Self-contained battery

Test Range

4-13 %A1C

Methodology

The A1CNow⁺ analyzer utilizes both immunoassay and chemistry technology to measure A1C.


CPT Codes

Tests performed with a CardioChek Plus professional analyzer and/or the A1CNow⁺ system may be reimbursable through Medicare and Medicaid.

SERVICE	CPT CODE
Capillary Blood Draw	36416
Hemoglobin A1C (Glycosylated Hemoglobin Test)	83036 QW
Glucose (Glucose; quantitative, blood (except reagent strip)	82947 QW
Lipid Panel (Total Cholesterol, HDL Cholesterol, Triglycerides)	80061 QW

Source: https://www.cms.gov/license/ama?file=/files/zip/ 22CLABQ1.zip . Accessed Oct. 3, 2022.

Test Strips, Controls, and Accessories



Self-wicking technology

- Plastic casing helps decrease breakage
- Available in other sample volumes

Ask about PreVantage

ptsdiagnostics.com/prevantage



A network of healthcare partnerships that close gaps in care

PreVantage solutions focus on...

- Population health improvement
- Improved efficiency
- Healthcare cost containment
- Patient satisfaction

Connectivity and Printing

PTS Connect® Printer

· Outputs results on paper or self-adhesive labels

- Connects to CardioChek analyzers through a provided cable
- Connects to A1CNow⁺ through PTS Connect Blu-Dock
- Includes one roll of self-adhesive labels



REF 2750



A1C Test Kits: 10-ct and 20-ct



A1CNow⁺ Professional Test System (10-cttestkit) REF 3024 A1CNow⁺ Professional Test System (20-ct test kit) REF 3021

Where to Buy

Through its large distribution network, PTS Diagnostics serves patients and healthcare providers all over the world.

For more information, schedule a call with one of our sales representatives at ptsdiagnostics.com/where-to-buy



CEas